

## **MDLand iClinic EHI Export**

### **§170.315(b)(10) Electronic Health Information Export**

MDLand iClinic is now in compliance with §170.315(b)(10) Electronic Health Information Export Certification requirements which specify the following:

- Enable a user to timely create an export file(s) with all of a single patient's electronic health information stored at the time of certification by the product, of which the Health IT Module is a part.
- A user must be able to execute this capability at any time the user chooses and without subsequent developer assistance to operate.
- Limit the ability of users who can create export file(s) in at least one of these two ways: (1) To a specific set of identified users (2) As a system administrative function, patient population electronic health information export- Create an export of all the electronic health information that can be stored at the time of certification by the product, of which the Health IT Module is a part.
- The export files(s) created must be electronic and in a computable format.
- The publicly accessible hyperlink of the export's format must be included with the exported file(s).

### **Export File Format in iClinic**

Patient EHI (both single patient and patient population) are exported in a zip folder in XML format. Please see XML format standard in <https://www.w3.org/TR/xml/>

### **A Detailed Guide on Exporting XML**

#### User privileges

Unauthorized users will not have access to "Advanced Settings" from where we can export the EHI. This is an existing functionality and is in line with § 170.315(b)(10) criteria.

#### Single Patient & Patient Population EHI Export

Authorized users can perform EHI export for either single patient or a specific patient population.

Step 1: Log in to the Account as a user with export privilege.

Step 2: Navigate to Settings/Advanced/ and Select 'Patient Data Export'. A user can select one of the options for EHI export and click "Submit".

- All Patients
- All Active Patients
- Patient with Encounter
- Patient with DOB
- Specify Patient [Single Patient]

**Settings**

- General
- EMR
- VBP Setting
- Patient Portal
- Billing/Account
- Schedule
- Follow Up
- ePrescription
- Lab
- Other Settings
- Advanced (ADM Only)
  - Privilege
  - Security Log
  - Audit Report
  - User List
  - My Account
  - MD Network
  - Patient Data Export**
  - Capitation Provider Setup
  - E-Payment

### Patient Data Export

Refresh New Close

Request Status: All

**Batch Editor**

Batch

☐ All Patients (Include Active, Inactive and Deceased)  
☐ All Active Patients  
☐ Patient with Encounter  
☐ Patient with DOB  
☐ Specify Patient

Submit Close

Step 3: After MDLand IT completes the processing of Patient EHI data file generation (it takes 10 minutes to a few hours depending on the data volume), the user can see a EHI export record for downloading and receive an iClinic message in his Inbox containing unzip password.

**Patient Data Export**

Refresh New Close

Request Status: All

Batch ID	Create Date	Status	Type	Patient List
1007	11/30/2023	Downloaded	Specify Patient	<a href="#">Detail</a>
1006	07/27/2023	Finished	All Active Patients	<a href="#">Detail</a>   <a href="#">Download</a>
1001	11/12/2020	Finished	All Active Patients	<a href="#">Detail</a>

Patient data is packed in a encrypted zip file. Password will be sent to you inbox when you download it.  
One task can be download once, please save the zip file in a safe place.

Step 4: Click “Download” to download the zip file and use the password received in the Inbox to unzip it and get the designation XML file.